2R2, 2R3, 2R37, 2R38, 4R21, 4R21=1, 4R52, 4R52=1

DE	Gebrauchsanweisung	4
EN	Instructions for use	11
FR	Instructions d'utilisation	19
IT	Istruzioni per l'uso	28
ES	Instrucciones de uso	36
PT	Manual de utilização	44
NL	Gebruiksaanwijzing	52
SV	Bruksanvisning	60
DA	Brugsanvisning	68
NO	Bruksanvisning	75
FI	Käyttöohje	83
PL	Instrukcja użytkowania	90
HU	Használati utasítás	98
CS	Návod k použití	106
RO	Instrucțiuni de utilizare	114
HR	Upute za uporabu	122
SL	Navodila za uporabo	130
SK	Návod na používanie	137
BG	Инструкция за употреба	145
TR	Kullanma talimatı	153
EL	Οδηγίες χρήσης	161
RU	Руководство по применению	169
JA	取扱説明書	178
ZH	使用说明书	185
KO	사용 설명서	192

1

4R21
4R21=1

4R52
4R52=1

2R37
2R38


2R2
2R3




2

3

INFORMATION

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INFORMATION

Date of the last update: 2014-10-24

- ▶ Please read this document carefully.
- ▶ Follow the safety instructions.

1.1 Construction and Function

The 2R2, 2R3, 2R37, 2R38 Tube Adapters and the 4R21, 4R21=1, 4R52 and 4R52=1 Tube Clamp Adapters are used as components for modular lower limb prostheses. The adapter combinations allow for controlled angle and length adaptation in the sagittal and frontal plane as well as adjustment of inward and outward rotation.

1.2 Combination Possibilities

For combination possibilities please consult the manufacturer.

2 Application

2.1 Intended Use

The product is intended solely for lower limb prosthetic fittings.

2.2 Area of Application

Area of application according to the MOBIS® mobility system:

2R2, 2R37, 4R21 and 4R52

Approved for a body weight of **up to 100 kg (220 lbs)**.

2R3, 2R38, 4R21=1 and 4R52=1

Approved for a body weight of **up to 136 kg (300 lbs)**.

The 4R21=1 and 4R52=1 Tube Clamp Adapters may only be used for TF prostheses and may only be positioned directly under the prosthetic knee joint or the prosthetic socket.

2.3 Environmental Conditions

Allowable Environmental Conditions
Temperature range for use: -10 °C to +60 °C (14 °F to 140 °F)
Allowable relative humidity 0 % to 90 %, non-condensing
Unallowable environmental conditions
Mechanical vibrations or impacts
Perspiration, urine, fresh water, salt water, acids
Dust, sand, highly hygroscopic particles (e. g. talcum)

2.4 Service Life

In general, all modular adapters are tested by the manufacturer for 3 million load cycles. Depending on the user's level of activity, this corresponds to a service life of 3 to 5 years.

3 Safety

3.1 Explanation of warning symbols

CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

3.2 General Safety Instructions

CAUTION

Using the product without following the instructions for use

Deteriorating health condition and damage to the product because of failure to follow the safety information

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with "**Inform the patient**".

CAUTION

Excessive strain on the product

Fall due to breakage of load-bearing components

- ▶ Use the prosthesis components according to the MOBIS classification (see section "Area of Application").
- ▶ **Inform the patient.**

CAUTION

Unallowable combination of prosthesis components

Fall due to breakage or deformation of the product

- ▶ Only combine the product with prosthesis components that are approved according to the section "Combination Possibilities".
- ▶ Based on the instructions for use of the prosthesis components, verify that they may be combined with each other.

⚠ CAUTION

Use under unallowable environmental conditions

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see section "Environmental Conditions").
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if applicable (e. g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).
- ▶ **Inform the patient.**

⚠ CAUTION

Exceeding the service life

Fall due to change in or loss of functionality and damage to the product

- ▶ Ensure that the approved service life is not exceeded (see section "Service Life").
- ▶ **Inform the patient.**

⚠ CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check a damaged product for proper functionality and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take suitable measures if required (e. g. repair, replacement, inspection by the manufacturer's customer service etc.).
- ▶ **Inform the patient.**

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthesis components relative to each other and by the development of noises.

4 Scope of Delivery

The scope of delivery for the product is shown on page 2 (Figure 1).

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

2R2, 2R3, 2R37, 2R38 Tube Adapters					
Fig.	Item no.		Quantity	Designation	Reference number
-	-	■	1	Instructions for Use	647H90=1
5	1	▲	4	Set screw	506G3=M8x12-V

4R21, 4R21=1, 4R52, 4R52=1 Tube Clamp Adapters					
Fig.	Item no.		Quantity	Designation	Reference number
-	-	■	1	Instructions for Use	647H90
4	1	▲	4	Set screw	506G3=M8x12-V
4	-	●	-	Single-components pack	4D4
	Consisting of:				
	2	-	1	Cap screw	-
	3	-	1	Straight pin	-
	4	-	1	Straight pin	-

For 4R21=1 and 4R52=1 Tube Clamp Adapters

4	5	■	1	Plastic ring	4X28=4
---	---	---	---	--------------	--------

For 4R21 and 4R21=1 Tube Clamp Adapters

4	6	■	1	Plastic ring	4X15=N2
---	---	---	---	--------------	---------

For 4R52 and 4R52=1 Tube Clamp Adapters

4	6	■	1	Plastic ring	4X15=N2-7
---	---	---	---	--------------	-----------

5 Preparation for Use

⚠ CAUTION

Incorrect alignment or assembly

Injuries due to fall resulting from damaged prosthesis components

- Observe the alignment and assembly instructions.

5.1 Adapting and Mounting the Tube Adapter

⚠ CAUTION

Incorrect processing of tube

Fall due to damage to the tube

- ▶ Do not clamp the tube into a vice.
- ▶ For shortening the tube, use only a tube cutter.

⚠ CAUTION

Incorrect mounting of the tube

Fall due to breakage of load-bearing components

- ▶ Slide the tube all the way to the stop in the intended prosthesis component when mounting.

⚠ CAUTION

Improper assembly of the screw connections

Fall due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the installation torques provided.
- ▶ Observe the instructions for securing the screw connections and the use of the correct length.

> **Recommended tools and materials:**

710D4 Torque Wrench, 719R3 Tube Cutter, 718R1 Tube Deburrer, degreasing cleaning agent (e.g. 634A3 Acetone)

- 1) Use the tube cutter to shorten the tube to fit the patient's measurements (see Fig. 2).
- 2) Carefully deburr the inside and outside of the tube with the deburring knife (see Fig. 3).
- 3) Loosen the cap screw (see Fig. 4, item 2) by approx. **2 turns**.
- 4) For definitive mounting, clean the insertion area of the connecting component on the inside and the tube on the outside using a degreasing cleaning agent.
- 5) Slide the tube **5 mm** into the connecting component while turning it slightly.
- 6) Take the connecting component with the inserted tube and place it on the work desk in a vertical position.
- 7) Slide the tube fully into the connecting component all the way to the stop.
- 8) Using the torque wrench, tighten the cap screw (for definitive mounting: **5 Nm**).

5.2 Mounting in the Modular Prosthesis

CAUTION

Improper assembly of the screw connections

Fall due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the installation torques provided.
- ▶ Observe the instructions for securing the screw connections and the use of the correct length.

INFORMATION

- ▶ Align the clamping slit in the tube clamp adapter to the anterior in order to avoid damage.

Mounting

> **Recommended tools and materials:**

710D4 Torque Wrench, 636K13 Loctite®

- 1) When installing, position the tube adapter distally and the connecting component proximally. For tube clamp adapters position the clamp slit anteriorly.
- 2) **For definitive mounting:**
Use Loctite® to secure the set screws (see Fig. 4 + 5, item 1).
- 3) Turn in the set screws.
- 4) **For trial fitting:**
Use the torque wrench to tighten the set screws (**10 Nm**).
- For definitive mounting:**
Pre-tighten the set screws with the torque wrench (**10 Nm**) and then tighten (**15 Nm**).
- 5) Use suitable set screws to replace any set screws that extend out too far or screw in too deeply (see selection table).

Selection Table for Set Screws

Reference number	Length (mm)
506G3=M8x12-V	12
506G3=M8x14	14
506G3=M8x16	16

Adjustment

The 4 set screws can be used to make static adjustments at any time during alignment, trial fittings, and even after the prosthesis is finished.

Replacement and Disassembly

When replacing or disassembling the prosthesis component, the previously adjusted position can be maintained if the two deepest screwed-in set screws, which are NOT located opposite to each other, are removed.

6 Maintenance

- The prosthesis components should be inspected after the first 30 days of use.
- Inspect the entire prosthesis for wear during normal consultations.
- Conduct annual safety inspections.

7 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return and collection processes.

8 Legal Information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

8.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

9 Technical Data

Reference number	2R2	2R3	4R21	4R21=1
Weight [g]	195	315	125	125
System height [mm]	–	–	33	33
Min. system height [mm]	97	97	–	–
Max. system height [mm]	232	472	–	–
Material	INOX (stainless steel)			
Diameter [mm]	30			
Max. body weight [kg]	100	136	100	136

Reference number	2R37	2R38	4R52	4R52=1
Weight [g]	160	275	80	80
System height [mm]	–	–	33	33
Min. system height [mm]	97	97	–	–
Max. system height [mm]	232	472	–	–
Material	Titanium			
Diameter [mm]	30			
Max. body weight [kg]	100	136	100	136

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Otto Bock HealthCare GmbH