## 12K41, 12K42, 12K44, 12K50

<b>DE</b> Gebrauchsanweisung (Benutzer) .....	2
<b>EN</b> Instructions for use (user) .....	13
<b>FR</b> Instructions d'utilisation (Utilisateur) .....	23
<b>IT</b> Istruzioni per l'uso (Utilizzatore) .....	34
<b>ES</b> Instrucciones de uso (Usuario) .....	45
<b>PT</b> Manual de utilização (Usuário) .....	56
<b>NL</b> Gebruiksaanwijzing (Gebruiker) .....	67
<b>SV</b> Bruksanvisning (Användare) .....	78
<b>RU</b> Руководство по применению (Пользователь) .....	89
<b>JA</b> 取扱説明書(ユーザー用) .....	101

## INFORMATION

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**INFORMATION**

Date of last update: 2020-10-08

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

The "12K41=\* ErgoArm, 12K42=\* ErgoArm plus, 12K44=\* ErgoArm Hybrid plus, 12K50=\* ErgoArm Electronic plus" product is referred to simply as the product/elbow joint below.


These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

**2 Product description**

**2.1 Construction and Function**

**2.1.1 Design**

1. Easy Plug
2. Strap clamp
3. Bowden cable (mechanical lock)
4. Dial
5. Humeral rotation feature adjusting screw

**2.1.2 Function**

The product is a mechanically/electronically driven elbow joint. It helps the wearer perform everyday tasks in combination with other prosthetic components (Combination possibilities). Various control programmes make optimal adaptation to the individual needs and abilities possible.

The ErgoArm models differ in the following characteristics:

Models	AFB	Easy Plug	Lock	Slip-stop
12K41=* ErgoArm	-	-	Mechanical	Mechanical

Models	AFB	Easy Plug	Lock	Slip-stop
<b>12K42=* ErgoArm plus</b>	x	-	Mechanical	Mechanical
<b>12K44=* ErgoArm Hybrid plus</b>	x	x	Mechanical	Mechanical
<b>12K50=* ErgoArm Electronic plus</b>	x	x	Mechanical or Electronic	Mechanical or Electronic

### 2.1.2.1 Definition of terms

#### AFB (Automatic Forearm Balance)

The AFB flexion aid stores the energy that is released when the prosthetic arm is extended and uses it to assist with flexion. The complex mechanism of the flexion aid controls the dynamics of forearm movement and equalises the leverage forces that develop via a compensation gear. As a result, the patient only has to exert minimal force on the flexion cable or through body momentum. The compensation coefficient can be easily adjusted to the individual weight of the prosthetic forearm and different clothing.

#### Lock

The product is either locked mechanically or electronically (see page 13). The internal ratchetless lock can be unlocked or locked in any position, even under a load. The product has a load capacity of up to 230 N in the locked state (this value may decline over the course of several years due to wear and tear). The lock will slip under higher loads.

#### Mechanical lock

A strong pull on the cable (audible clicking sound) locks or unlocks the joint.

Thanks to the integrated compensation mechanism, the lock can be released with little effort, even under load. The integrated strain relief reduces the risk of damage when the pull cable is pulled too hard.

#### Electronic lock

The lock can be released and engaged using electrode signals or switch signals (depending on the chosen switching version).

#### Slip-stop function

The slip-stop function makes it possible to lower the forearm in a controlled manner without having to fully release the lock.

#### Slip-stop function, mechanical

A gentle pull on the cable (no audible click) unlocks the product so the forearm can be moved continuously. When the pull on the cable stops, the product is locked again.

#### Slip-stop function, electronic

The slip-stop function is activated by the electrode signals (depending on the chosen control version). When the electrode signal stops, the product is locked again.

## 3 Intended use

### 3.1 Indications for use

The product is intended **exclusively** for upper limb exoprosthetic fittings.

### 3.2 Conditions of use

The product is designed for use by unilateral or bilateral amputees.

The product is intended **exclusively** for use on adults.

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, parachuting, paragliding, etc.).

### 3.3 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

### 3.4 Qualification

The fitting of a user with the product may only be carried out by O&P professionals who have been authorised through the corresponding Ottobock training.

## 4 Safety

### 4.1 Explanation of warning symbols

 <b>WARNING</b>	Warning regarding possible serious risks of accident or injury.
 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

### 4.2 Structure of the safety instructions

 <b>CAUTION</b>
<p><b>The heading describes the source and/or the type of hazard</b></p> <p>The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:</p> <ul style="list-style-type: none"><li>&gt; E.g.: Consequence 1 in case of failure to observe the hazard</li><li>&gt; E.g.: Consequence 2 in case of failure to observe the hazard</li></ul> <p>▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.</p>

### 4.3 General safety instructions

 <b>WARNING</b>
<p><b>Non-observance of safety notices</b></p> <p>Personal injury/damage to the product due to using the product in certain situations.</p> <p>▶ Observe the safety notices and the stated precautions in this accompanying document.</p>

 <b>WARNING</b>
<p><b>Operating the product near active implanted systems</b></p> <p>Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic interference of the product.</p> <ul style="list-style-type: none"><li>▶ When operating the product in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.</li><li>▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.</li></ul>

 <b>WARNING</b>
<p><b>Use of the product while operating a vehicle or machinery</b></p> <ul style="list-style-type: none"><li>&gt; Risk of accident due to unexpected behaviour of the product.</li><li>&gt; Risk of injury due to faulty control or malfunction of the product.</li></ul> <ul style="list-style-type: none"><li>▶ Applicable national legal regulations must be complied with when operating vehicles of any kind or when operating machinery with the product.</li><li>▶ Have your driving ability examined and verified by an authorised test centre.</li></ul>

- ▶ Switch off the product before driving.

**⚠ WARNING**

**Charging the prosthesis without taking it off**

Risk of electric shock due to defects in the power supply unit or in the battery charger.

- ▶ For safety reasons, remove the prosthesis prior to charging it.

**⚠ CAUTION**

**Changing prosthetic components when switched on**

Risk of injury due to faulty control or malfunction of the prosthesis solution.

- ▶ Remove the rechargeable battery from the battery receptacle or switch off the prosthesis solution before changing prosthetic components (such as the terminal device).
- ▶ The 12K50=\* ErgoArm can be locked and unlocked by means of the mechanical pull cable even when it is switched off.

**⚠ CAUTION**

**Overloading due to unusual activities**

Injury due to unexpected product behaviour as a result of a malfunction.

- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures your personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately by an O&P professional. If necessary, he or she will forward the product to an authorised Ottobock Service Centre.

**⚠ CAUTION**

**Incorrect donning and doffing of the prosthesis**

Risk of injury due to sudden flexion of the prosthesis.

- ▶ Switch off the prosthesis before donning or doffing.
- ▶ The prosthesis must always be flexed for donning and doffing.

**⚠ CAUTION**

**Sudden flexion of the product due to high level of compensation**

Risk of injury due to sudden flexion of the prosthesis.

- ▶ Note the level of compensation when taking off heavy clothing.
- ▶ Note the setting for the level of compensation when the prosthesis is not being worn.
- ▶ The prosthesis must always be flexed for donning and doffing.

**⚠ CAUTION**

**Manual unlocking of elbow lock under load**

Injury by release of elbow lock under load.

- ▶ Particular caution should be exercised when unlocking the elbow lock while lifting heavy loads.
- ▶ Be careful when unlocking the lock under such conditions due to the possibility of injury.

**⚠ CAUTION**

**Penetration of dirt and humidity into the product**

Risk of injury due to unexpected product behaviour or malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.
- ▶ Do not expose the product, and especially the elbow joint, to dripping or splashing water.
- ▶ Wear the product, and especially the elbow joint, under suitably resistant clothing in the rain.

**⚠ CAUTION**

**Signs of wear and tear on the product components**

Injury due to faulty control or malfunction of the product

- ▶ In the interest of preventing injuries and to maintain product quality, we recommend performing regular service.
- ▶ Please contact your O&P professional for more information.

**⚠ CAUTION**

**Unintentional unlocking of the terminal device**

Risk of injury due to releasing the terminal device from the forearm (e.g. while carrying objects).

- ▶ When using a quick-disconnect wrist unit, position the terminal device so that slight twisting does not release the terminal device from the forearm.

**⚠ CAUTION**

**Mechanical stress on the product**

Injury due to faulty control or malfunction of the product components.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

**⚠ CAUTION**

**Independent manipulation of the product**

Injury due to product damage or malfunction.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

**⚠ CAUTION**

**Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)**

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the product.

### **⚠ CAUTION**

#### **Distance to HF communication devices is too small (e.g. mobile phones, Bluetooth devices, WiFi devices)**

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping the following minimum distances to these HF communication devices is recommended:
  - Mobile phone GSM 850/GSM 900: 0.99 m
  - Mobile phone GSM 1800/GSM 1900/UMTS: 0.7 m
  - DECT cordless phones incl. base station: 0.35 m
  - WiFi (routers, access points,...): 0.22 m
  - Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

### **⚠ CAUTION**

#### **Risk of pinching in the joint flexion area**

Injuries due to pinching of body parts.

- ▶ Ensure that fingers and other body parts are not in this area when bending the joint.

### **⚠ CAUTION**

#### **Charging the product with damaged power supply unit/charger/charger cable**

Injury due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Check the power supply unit, charger and charger cable for damage before use.
- ▶ Replace any damaged power supply unit, charger or charger cable.

### **NOTICE**

#### **Coating, gluing or painting the prosthesis**

Damage or fracture due to chemical processes.

- ▶ The prosthesis must not be coated, glued or painted.

### **NOTICE**

#### **Improper product care**

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Clean the product and the contacts with a damp cloth and mild soap only.

### **NOTICE**

#### **Damage to the battery**

Damage to the battery due to using both connection options.

- ▶ Only use one of the two connections for the battery (connection on the forearm cable or Easy Plug connection).

## **5 Scope of Delivery and Accessories**

### **12K41 ErgoArm plus scope of delivery**

- 1 pc. 12K41=\* ErgoArm
- 1 pc. 647G1475 Instructions for use (user)

### **12K42 ErgoArm plus scope of delivery**

- 1 pc. 12K42=\* ErgoArm plus
- 1 pc. 647G1475 Instructions for use (user)

### **12K44 ErgoArm Hybrid plus scope of delivery**

- 1 pc. 12K44=\* ErgoArm Hybrid plus
- 1 pc. 647G1475 Instructions for use (user)

## 12K50 ErgoArm Electronic plus scope of delivery

- 1 pc. 12K50=\* ErgoArm Electronic plus
- 1 pc. 647G1475 Instructions for use (user)

## 6 Use

### 6.1 Adjusting the level of compensation

Complete the following steps to set the level of compensation:


- 1) Flex the arm.
- 2) Use the dial to adjust the compensation.

**INFORMATION: The adjustment range is limited by stops.**

**INFORMATION: The adjustment mechanism is protected by a slip clutch.**

There are three different compensation ranges:

#### **Range 1 = low compensation level**

Permits the arm to swing freely while walking.

#### **Range 2 = moderate compensation level**

Balance is positively supported when flexing the arm. Compensation is reduced again during extension.

#### **Range 3 = constant compensation level**

If set correctly, the weight of the forearm is balanced by the compensation level.

#### **Range 4 = low compensation level**

Low compensation until flexion stop.

### 6.2 Adjusting the upper arm rotation joint

Complete the following steps to rotate and adjust the upper arm rotation joint:

- 1) Loosen the adjusting screw.
- 2) Move the upper arm rotation joint to the desired position.
- 3) Re-tighten the adjusting screw.

## 7 Charging the battery

The following points must be observed when charging the battery:

- The capacity of a fully charged battery is sufficient for one full day.
- We recommend charging the product once a day when used by the patient on a daily basis.
- When used daily, the complete charging unit (power supply – battery charger) may remain plugged into the wall socket.

### 7.1 757L20 battery charger and 757B2\* EnergyPack

#### NOTICE

#### Deep discharge of the battery

Irreparable damage to the product and resulting loss of functionality in case of irregular and incomplete charging of the battery.

- ▶ In case of long-term storage, fully charge the battery every 4 to 6 months.

- 1) Plug the charging station with power supply unit into the wall socket.  
→ The green LED on the power supply unit and the green LED in the middle of the LED row light up.
- 2) Insert one or two rechargeable batteries into the charging station.
- 3) The green LED for the respective charging socket lights up ● and the charging process starts.
- 4) When the rechargeable battery is fully charged, the green LED for the respective charging socket flashes ✨.
- 5) Remove the rechargeable battery after the charging process is complete.

### 7.2 757L35 battery charger and 757B35=\* rechargeable battery

- 1) Connect the charging plug to the charging receptacle on the prosthetic socket.  
→ A correct connection between the battery charger and the prosthesis is indicated by feedback (Status signals).
- 2) The charging process starts.  
→ The prosthesis is switched off automatically.
- 3) Disconnect the prosthesis when the charging process is complete.

## 8 Cleaning

- 1) Clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock Derma Clean) when needed.  
Make sure that no liquids get into the product and product components.
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

## 9 Maintenance

Performing regular maintenance (service inspections) every 24 months is recommended to prevent injuries and maintain the quality of the product.

The grace period is no more than one month before or three months after maintenance is due.

In general, all products are subject to compliance with the maintenance intervals during the warranty period. This is the only way to maintain full warranty cover.

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

The following components always have to be submitted to the O&P professional for maintenance and repairs:

The prosthesis, battery charger and power supply unit.

## 10 Legal information

### 10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 10.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

### 10.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

## 11 Technical data

Ambient conditions	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing
Charging the battery	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing

Reference number	12K41=* 12K42=* 12K44=* 12K50=*	12K41=* 12K42=* 12K44=* 12K50=*	12K41=* 12K42=* 12K44=* 12K50=*	12K41=* 12K42=* 12K44=* 12K50=*
Current draw (in mA)	-	-	-	3.5
Power supply	-	-	-	757B20 757B21 757B25=* 757B35=*
Max. load in locked state (N)	230			
Flexion angle (in °)	15°-145°			
Service life	5 years			

## 12 Appendices

### 12.1 Symbols Used

Declaration of conformity according to the applicable European directives

Serial number (YYYY WW NNN)  
YYYY – year of manufacture  
WW – week of manufacture  
NNN – sequential number

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Manufacturer

Medical device

### 12.2 Display of the current charge level

#### 12.2.1 757B2\* EnergyPack

LED indicator	Event
	Charge level above 50% (illuminated in green)
	Charge level below 50% (flashing alternating green and orange)
	Charge level below 5% (illuminated in orange)

When inserting/engaging the rechargeable battery in the prosthesis, a charge level indicator is activated on the battery for a few seconds.

#### 12.2.2 757B35=\* Battery

The battery charge level can be queried at any time.

- 1) With the prosthesis switched on, press the charging receptacle button and hold for less than one second.
- 2) The LED display on the charging receptacle provides information on the current battery charge level (Status signals).

Charging receptacle	Event
	Battery fully charged (illuminated in green)
	Battery 50% charged (illuminated in yellow)
	Battery drained (illuminated in orange)

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