



## 17PK1=\*


<b>DE</b>	Gebrauchsanweisung .....	4
<b>EN</b>	Instructions for use .....	10
<b>FR</b>	Instructions d'utilisation .....	16
<b>IT</b>	Istruzioni per l'uso .....	23
<b>ES</b>	Instrucciones de uso .....	29
<b>PT</b>	Manual de utilização .....	36
<b>NL</b>	Gebruiksaanwijzing .....	42
<b>SV</b>	Bruksanvisning .....	49
<b>DA</b>	Brugsanvisning .....	54
<b>NO</b>	Bruksanvisning .....	60
<b>PL</b>	Instrukcja użytkowania .....	66
<b>HU</b>	Használati utasítás .....	72
<b>CS</b>	Návod k použití .....	79
<b>HR</b>	Upute za uporabu .....	84
<b>TR</b>	Kullanma talimatı .....	90
<b>RU</b>	Руководство по применению .....	97

2

743L100

## INFORMATION

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### INFORMATION

Last update: 2017-03-15

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These Instructions for Use provide you with important information on the installation and processing of the 17PK1=\* System Knee Joints made of plastic material.

## 2 Product description

### 2.1 Components/design

Item	Component	Reference number			
		17PK1=14		17PK1=20	
1	Countersunk head screw	29PK1=14		29PK1	
2	Lock cover				
3	Unlock cable (perlon)				
4	Compression spring				

Item	Component	Reference number			
		17PK1=14		17PK1=20	
5	Threaded sleeve	29PK1=14	21A12	29PK1	21A12
6	Lock				
7	Upper joint section				
8	Receiver for set screw				
9	M4x12 set screw	506G21=M4x12		506G21=M4x12	
10	Bearing bushing	29PK2=14		29PK2	
11	Joint bolt				
12	Countersunk head screw (17PK1=14) Cap screw (17PK1=20)				
13	Joint nut with M6 internal thread				
14	Lower joint section				
15	System bar receiver				
16	Countersunk head screw (17PK1=14) Cap screw (17PK1=20)	29PK4=14		29PK4	
17	Parallel truing tube adapter (17PK1=14 only)				
18	Parallel truing tube				
19	System bar (not included in the scope of delivery!)				

### 3 Intended use

#### 3.1 Indications for use

The system knee joints are to be used **exclusively** in pairs for the orthotic fitting of the lower limbs on a patient.

Reference number	Max. body weight
17PK1=14	45 kg
17PK1=20	100 kg

#### 3.2 Indications

Partial or total paralysis of the leg muscles.  
Indications must be determined by the physician.

#### 3.3 Further Usage Restrictions

The product is not suitable for use in sports.

## 4 Safety

### 4.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 4.2 Safety instructions

 **CAUTION**

**Excessive strain on load-bearing components**

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service etc.).

 **CAUTION**

**Incorrect alignment or assembly**

Injuries due to changes in or loss of functionality

- ▶ The product may only be installed by trained, qualified personnel.
- ▶ Observe the alignment and assembly instructions.

 **CAUTION**

**Mechanical damage to the product**

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

 **CAUTION**

**Excessive strain due to use on more than one patient**

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

**⚠ CAUTION****Unallowable use of lubricants**

Risk of injury and loss of functionality as well as damage to the product

- ▶ The product is designed to be free of lubricants. Avoid contact with lubricants (e.g. grease, oil).
- ▶ After contact with lubricants, promptly clean the product or parts of the product with a degreasing cleaning agent (e.g. acetone or isopropyl alcohol).

**NOTICE****Contact with saltwater or water containing chlorine/soap, fresh water or abrasive substances (e.g. sand)**

Damage and premature product wear

- ▶ Following any contact with the substances identified above, promptly clean the product in accordance with the section "Cleaning".

**NOTICE****Exposure of the product to unsuitable environmental conditions**

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

## 5 Preparation for use

**INFORMATION**

- ▶ Check the exact parallel alignment of the system joints before delivering the orthosis.

**INFORMATION**

In order to facilitate the parallel alignment of the system joints, each pair of system joints is supplied with an adjustment aid. The adjustment aid consists of 2 screws, 1 parallel truing tube and, depending on the size of the system joints, 2 adapter pieces.

Due to varying alignment parameters, such as the generated forefoot leverage or the flexion position of the orthosis, the system knee joints are subjected to increased strain.

### **Assembling the lock**

The lock is pre-assembled in the orthosis joint upon delivery. The 29PK1\* set is required to replace the lock. Insertion aids are attached to the locks so that trying them out is easier.

- 1) Insert the bearing bushings (10) into the upper joint section (7).
- 2) Slide the lower joint section (14) into the upper joint section (7) and secure it with the joint bolt (11), the joint nut with internal thread (13) and the cap screw (12).
- 3) Insert the lock (6) in size 3 into the lock channel of the joint up to the stop.
- 4) Check the play of the joint.

**INFORMATION: If the joint has play or the lock cannot be inserted up to the stop, exchange the lock for a different size (1, 2, 4 or 5).**

- 5) Break the chosen lock off the insertion aid.
- 6) Screw the threaded sleeve (5) onto the unlock cable (3).
- 7) Insert the unlock cable (3) with threaded sleeve (5) into the chosen lock (6).
- 8) Guide the unlock cable (3) through the compression spring (4) and lock cover (2).
- 9) Insert the lock (6) and the compression spring (4) with the unlock cable (3) into the lock channel.
- 10) Attach the lock cover (2) in the upper joint section (7) with the counter-sunk head screw (1).

### **Installing the system bar in the system knee joint**

- 1) Insert the system bar (19) into the system bar receiver (15) and mark the bore on the system bar (19).
- 2) Remove the system bar (19) from the system bar adapter (15) and drill a **Ø 4 mm** bore in the system bar (19).
- 3) Secure the system bar (19) and the system knee joint with the set screw (9).

### **Gluing the system bar**

For stability reasons, the system bar (19) must be glued into the system bar receiver (15) of the system knee joints using 636W28 Ottobock special adhesive.

- 1) Use a degreasing cleaner (such as 634A1) to degrease the bonding surfaces of the system bar (19) and the system bar receiver (15).

- 2) Mix the special adhesive and hardener thoroughly and apply it to the components being glued.
- 3) Insert the system bar (19) and secure it with the set screw (9). The final bonding strength of the adhesive is reached after **16 h**.
- 4) Prior to delivering the orthosis:  
Secure the joint bolt (11) and the joint nut (13) with 636K13 Loctite 241.  
Verify the proper functionality of the product.

### **Removing the system bar**

- 1) Unscrew the set screw (9).
- 2) Tap out and remove the system bar (19).

## **6 Cleaning**

After contact with water containing salt, chlorine or soap, or when they get dirty, the system joints must be cleaned promptly:

- 1) Rinse the system joints with clean fresh water.
- 2) Dry the system joints with a cloth or allow to air dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

The system joints must be cleaned promptly after coming into contact with lubricants (e.g. grease, oil):

- 1) Moisten a cloth with a degreasing cleaning agent (e.g acetone or isopropyl alcohol).
- 2) Wipe down the system joints or components with the cloth.

## **7 Maintenance**

### **⚠ CAUTION**

#### **Use of an incorrect lock**

Injuries from falling due to an incorrect lock being installed

- ▶ Only use locks of the same size when reinstalling or replacing.

### **ℹ INFORMATION**

The lock is chosen according to the joint and may not be used for installation on other joints.

- ▶ Observe the information in the service set when replacing parts.

The manufacturer recommends a semi-annual inspection to verify functionality and check for wear of the system joints.

Spare parts are listed in the section "Components/design".

## 8 Disposal

Dispose of the product in accordance with national regulations.

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 9.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH