

Evidence Essentials.

C-Leg microprocessor knee

	Mobility need or deficit of the patient	Evidence for benefits of the C-Leg vs. NMPK
Safety	<p>Patient has a history of fall-related injury</p> <p>Patient stumbles and/or falls repeatedly</p> <p>Patient avoids activities due to fear of falling</p> <p>Patient sustained fall-related injuries</p>	<ul style="list-style-type: none"> - Significant 65% reduction of patients who experience injurious falls (Campbell et al., 2020) - Significant reduction in falls of up to 80% (Hafner et al., 2007 and 2009; Highsmith et al., 2010; Kahle et al., 2008; Kannenberg et al., 2014, Kaufman et al., 2018; Davie-Smith et al., 2021) - Significant reduction in stumbles of up to 57% (Hafner et al., 2009; Highsmith et al., 2010; Kahle et al., 2008; Kannenberg et al., 2014) - Significant improvements in balance and indicators for the risk of falling, such as Timed-up-and-go-test, forced gait perturbations in the gait lab, ABC scale, etc. (Blumentritt et al., 2009; Burnfield et al., 2012; Hafner et al., 2007 and 2009; Kannenberg et al., 2014; Kaufman et al., 2007; Lansade et al., 2018; Davie-Smith et al., 2021)
Mobility	<p>Patient has difficulty negotiating slopes/hills</p>	<ul style="list-style-type: none"> - Significant improvement in quality of slope/hill descent towards natural, reciprocal (step-over-step) gait pattern (Hafner et al.; 2007 and 2009; Highsmith et al., 2013; Kannenberg et al., 2014) - Significant increase in downhill walking speed of up to 40% (Burnfield et al., 2012; Hafner et al., 2007 and 2009; Highsmith et al., 2013; Kannenberg et al., 2014)

	Mobility need or deficit of the patient	Evidence for benefits of the C-Leg vs. NMPK
Mobility	Patient has difficulty negotiating uneven terrain and obstacles	<ul style="list-style-type: none"> - Significant increase in walking speed on uneven terrain and obstacle courses of up to 21% (Kahle et al., 2008; Seymour et al., 2007)
Mobility	Patient has difficulty descending stairs with reciprocal (step-over-step) gait	<ul style="list-style-type: none"> - Significant improvement in quality of stair descent towards natural, reciprocal (step-over-step) gait pattern (Hafner et al., 2007 and 2009; Kahle et al., 2008; Kannenberg et al., 2014; Schmalz et al., 2007)
Mobility	Patient has difficulty with dual tasking while walking with the prosthesis	<ul style="list-style-type: none"> - Significant reduction in cognitive demand while walking with the prosthesis (Hafner et al., 2007 and 2009; Seymour et al., 2007; Williams et al., 2006) - Significant reduction in cortical brain activity and perfusion during dual-tasking (Möller et al., 2019; Ramstrand et al., 2020) - Significantly improved capacity and performance in executing a concurrent task while walking with the prosthesis (Hafner et al., 2007 and 2009; Morgan et al., 2015; Seymour et al., 2007; Williams et al., 2006)
Mobility	Patient is limited in his/her mobility	<ul style="list-style-type: none"> - About 50% of K2 patients are able to improve their overall mobility level to K3 (Hafner et al. 2009; Kahle et al., 2008; Kannenberg et al., 2014)
Musculo-skeletal pain	Patient suffers from joint and back pain due to gait asymmetry and excessive loading	<ul style="list-style-type: none"> - Significant improvement in gait symmetry and, thus, loading of the locomotor system (Kaufman et al, 2007 and 2012; Segal 2006) - Significant stance knee flexion that results in shock absorption to unload proximal joints and the spine (Kaufman et al., 2007; Segal et al., 2006)

References

- Blumentritt S, Schmalz T, Jarasch R. The safety of C-leg: Biomechanical tests. J Prosthet Orthot 2009;21(1):2-17. [Download](#)
- Burnfield JM, Eberly VJ, Gronely JK, Perry J, Yule WJ, Mulroy SJ. Impact of stance phase microprocessor-controlled knee prosthesis on ramp negotiation and community walking function in K2 level transfemoral amputees. Prosthet Orthot Int 2012;36(1):95-104. [Download](#)
- Campbell JH, Stevens PM, Wurdeman SR. OASIS I: Retrospective analysis of four different microprocessor knee types. Journal Rehabil Assist Technol Eng 2020;7: 1-10. [Download](#)
- Davie-Smith F, Carse B. Comparison of patient-reported and functional outcomes following transition from mechanical to microprocessor knee in the low-activity user with a unilateral transfemoral amputation. Prosth Orthot Int 2021;45(3):198-204. [Download](#)
- Hafner BJ, Willingham LL, Buell NC, Allyn KJ, Smith DG. Evaluation of Function, Performance, and Preference as Transfemoral Amputees Transition from Mechanical to Microprocessor Control of the Prosthetic Knee. Arch Phys Med Rehabil 2007;88(2):207-17. [Download](#)
- Hafner BJ, Smith DG. Differences in function and safety between Medicare Functional Classification Level-2 and -3 transfemoral amputees and influence of prosthetic knee joint control. J Rehabil Res Dev 2009;46(3):417-434. [Download](#)
- Highsmith MJ, Kahle JT, Bongiorno DR, Sutton BS, Groer S, Kaufman KR. Safety, energy efficiency, and cost efficacy of the C-leg for transfemoral amputees. Prosthet Orthot Int 2010;34(4):362-377. [Download](#)
- Highsmith MJ, Kahle JT, Miro RM, Mengelkoch, MJ: Ramp descent performance with the C-leg and interrater reliability of the Hill Assessment Index. Prosthet Orthot Int 2013; 37(5): 362-368. [Download](#)
- Kahle JT, Highsmith MJ, Hubbard SL. Comparison of Non-microprocessor Knee Mechanism versus C-Leg on Prosthesis Evaluation Questionnaire, Stumbles, Falls, Walking Tests, Stair Descent, and Knee Preference; J Rehabil Res Dev 2008;45(1):1-14. [Download](#)
- Kannenbergh A, Zacharias B, Pröbsting E: Benefits of microprocessor prosthetic knees to limited community ambulators: A systematic review. J Rehabil Res Dev 2014;51(10):1469-1495. [Download](#)
- Kaufman KR, Levine JA, Brey RH, et al. Gait and Balance of transfemoral amputees using passive mechanical and microprocessor-controlled prosthetic knees. Gait Posture 2007;26:489-493. [Download](#)
- Kaufman KR, Frittoli S, Frigo CA. Gait asymmetry of transfemoral amputees using mechanical and microprocessor controlled prosthetic knees. Clin Biomech 2012;27(5):460-465. [Download](#)

Kaufman KR, Bernhardt KA, Symms K. Functional assessment and satisfaction of transfemoral amputees with mobility (FASTK2): A clinical trial of microprocessor-controlled vs. non-microprocessor-controlled knees. Clin Biomech (Bristol, Avon) 2018 Oct;58:116-122. [Download](#)

Lansade C, Vicaut E, Paysant J, Ménager D, Cristina MC, Braatz F, Domayer S, Pérennou D, Chiesa G. Mobility and safety with a microprocessor-controlled knee in moderately active amputees: A multi-centric randomized crossover trial. Ann Phys Rehabil Med 2018;61(5):278-285. [Download](#)

Morgan SJ, Hafner BJ, Kelly VE. The effects of a concurrent task on walking in persons with transfemoral amputation compared to persons without limb loss. Prosthet Orthot Int 2016 Aug;40(4):490-496. [Download](#)

Möller S, Rusaw D, Hagberg K, Ramstrand N. Reduced cortical brain activity with the use of microprocessor-controlled prosthetic knees during walking. Prosthet Orthot Int 2019;43(3):257-265. [Download](#)

Ramstrand N, Rusaw DF, Möller SF. Transition to a microprocessor controlled prosthetic knee: Executive functioning during single and dual-task gait. Prosthet Orthot Int 2020;44(1):27-35. [Download](#)

Segal AD, Orendurff MS, Klute GK, McDowell ML, Pecoraro JA, Shofer J, Czerniecki JM. Kinematic and kinetic comparisons of transfemoral amputee gait using C-Leg and Mauch SNS prosthetic knees. J Rehabil Res Dev 2006;43(7):857-870. [Download](#)

Seymour R, Engbretson B, Kott K, Ordway N, Brooks G, Crannell J, Hickernell E, Wheller K. Comparison between the C-leg(R) microprocessor-controlled prosthetic knee and non-microprocessor control prosthetic knees: A preliminary study of energy expenditure, obstacle course performance, and quality of life survey. Prosthet Orthot Int 2007;31(1):51- 61. [Download](#)

Schmalz T, Blumentritt S, Marx B. Biomechanical Analysis of Stair Ambulation in Lower Limb Amputees. Gait Posture 2007;25:267-278. [Download](#)

Williams RM, Turner AP, Orendurff M, Segal AD, Klute GK, Pecoraro J, Czerniecki J. Does Having a Computerized Prosthetic Knee Influence Cognitive Performance during Amputee Walking? Arch Phys Med Rehabil 2006;87:989-994. [Download](#)

Ottobock Reimbursement North America

P 800 328 4058 F 800 230 3962

US: shop.ottobock.us

CA: shop.ottobock.ca

reimbursement911@ottobock.com