





50C91 Smartspine Universal Collar

DE	Gebrauchsanweisung	4
EN	Instructions for use	8
FR	Instructions d'utilisation	13
IT	Istruzioni per l'uso	18
ES	Instrucciones de uso	23
PT	Manual de utilização	28
NL	Gebruiksaanwijzing	32
SV	Bruksanvisning	37
DA	Brugsanvisning	42
NO	Bruksanvisning	46
FI	Käyttöohje	51
IS	Notkunarleiðbeiningar	55
PL	Instrukcja użytkowania	60
HU	Használati utasítás	64
CS	Návod k použití	69
RO	Instrucțiuni de utilizare	74
HR	Upute za uporabu	78
SK	Návod na používanie	83
BG	Инструкция за употреба	87
TR	Kullanma talimatı	92
JA	取扱説明書	97
ZH	使用说明书	101
KO	사용 설명서	105

	Größe / Size	Abstand* / Distance	
		cm	inch
50C91	universal	6,4–12,7	2,5–5

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INFORMATION

Last update: 2018-08-10

► Please read this document carefully before using the product.

- | |
|---|
| <ul style="list-style-type: none">▶ Follow the safety instructions to avoid injuries and damage to the product.▶ Instruct the user in the proper and safe use of the product.▶ Please keep this document in a safe place. |
|---|

These instructions for use provide important information on the fitting and application of the 50C91 Smartspine Universal Collar cervical brace.

2 Intended use

2.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the cervical spine and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

2.2 Indications

50C91 Smartspine Universal Collar

- Acute pain syndromes of the cervical spine
- Stenosis of the cervical spine
- Whiplash trauma of the cervical spine
- Fractures in the cervical spine (simple and stable)
- Dislocation & dislocation fractures of the cervical spine (postoperative)
- Traumatic injury to the cervical spine
- Radiculopathies of the cervical spine

The indication must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

None known.

2.3.2 Relative Contraindications

Consult the physician if any of the following occur: skin diseases/injuries, inflammation, raised scars with swelling, reddening, and warming in the area where the brace is fitted, tumors, lymphatic congestion and dysesthesia in the cervical spine.

2.4 Mechanism of Action

The brace immobilizes the cervical spine.

3 Safety

3.1 Explanation of warning symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

3.2 General safety instructions

CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

NOTICE

Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

NOTICE

Contact with ointments, lotions, or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to ointments, lotions, or other products that contain oils or acids.

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Selecting the Size

The brace is available in a universal size. Distances between the tip of the chin and the top of the sternum in the range of 6.4 cm to 12.7 cm can be fitted.

4.2 Adaptation

The brace is delivered ready for trial fitting. Adjustment is not required.

- 1) Turn the brace lock on both sides counter-clockwise to the "open" position (arrowheads in vertical position) (see fig. 1).
- 2) Engage the brace on both sides in the corresponding position (1-9) (see fig. 2).

CAUTION! Ensure that the brace setting is the same on both sides.

- 3) Turn the brace lock clockwise to the "closed" position (arrowheads in horizontal position) (see fig. 3).
- 4) **Adjusting the inner part of the brace for small neck diameters:** Hold the neckpiece of the inner part of the brace with both hands and press it together. This will allow you to bend the brace round.
- 5) **Optional:** Shorten the hook-and-loop closure to the desired length (see fig. 6).

4.3 Application

⚠ CAUTION

Use with heat-sensitive skin

Skin irritation due to overheating

- ▶ Do not use the product in case of a known heat allergy.
- ▶ Do not continue to use the product if skin irritation occurs.

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the brace is applied properly and fits correctly.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

- > The patient should be standing or sitting down.
- > Open the hook-and-loop closure on one side of the brace.
- 1) Place the inside of the brace on the neck from the front (see fig. 4).
- 2) Place the outside of the brace on the neck from behind.
 - The chin rests on the recess.
- 3) Fasten the hook-and-loop closure again (see fig. 5).

5 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the brace regularly:

- 1) Fasten all hook-and-loop closures.
- 2) Hand wash the brace in warm water at **30 °C (86 °F)** with standard mild detergent.
- 3) Rinse well.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

6 Disposal

Dispose of the product in accordance with national regulations.

7 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

7.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

7.2 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

7.3 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The manufacturer therefore drew up the declaration of conformity on its own responsibility in accordance with Annex VII of the directive.

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