



Evidence Essentials

Kenevo/Microprocessor Knees for K2

	Mobility need or deficit of the patient	Evidence for benefits of Kenevo/MPK vs. NMPK in K2 patients
Safety	<p>Patient stumbles and/or falls repeatedly</p> <p>Patient avoids activities due to fear of falling</p> <p>Patient sustained fall-related injuries</p>	<ul style="list-style-type: none"> - Significant reduction in falls of up to 80% (Hahn et al., 2021; Davie-Smith et al., 2021; Kaufman et al., 2018; Mileusnic et al., 2017; Wong et al., 2015; Hahn et al., 2015; Kannenberg et al., 2014; Hafner et al., 2009; Kahle et al., 2008) - Significant reduction in fear of falling (Hahn et al., 2021; Jayaraman et al., 2021; Mileusnic et al., 2017; Wong et al., 2015; Hahn et al., 2015) - Significant reduction in the frequency of stumbles (Mileusnic et al., 2017; Kannenberg et al., 2014; Hafner et al., 2009) - Significant improvements in balance and indicators for the risk of falling, such as Timed-up-and-go-test, ABC scale, PEQ Addendum; Modified Falls Efficacy Scale, etc. (Hahn et al., 2021; Davie-Smith et al., 2021; Jayaraman et al., 2021; Lansade et al., 2018; Hahn et al., 2016; Wong et al., 2015; Kannenberg et al., 2014; Burnfield et al., 2012; Hafner et al., 2007 and 2009)
Mobility	Patient has difficulty negotiating slopes/hills	<ul style="list-style-type: none"> - Significant improvement in quality of slope descent towards more natural gait pattern (Kannenberg et al., 2014; Burnfield et al., 2012; Hafner et al., 2009) - Significant increase in downhill walking speed of up to 36% (Kannenberg et al., 2014; Burnfield et al., 2012; Hafner et al., 2009) - Significant improvement in patient-reported slope ambulation (Hahn et al., 2016)
Mobility	Patient has difficulty negotiating uneven terrain and obstacles	<ul style="list-style-type: none"> - Significant increase in walking speed on uneven terrain and obstacle courses of up to 20% (Kannenberg et al., 2014; Hafner et al., 2009; Kahle et al., 2008) - Significant improvement in patient-reported uneven terrain and obstacle negotiation (Hahn et al., 2016)
Mobility	Patient has difficulty descending stairs with reciprocal (step-over-step) gait	<ul style="list-style-type: none"> - Significant improvement in quality of stair descent towards more natural gait pattern (Kannenberg et al., 2014; Hafner et al., 2009; Kahle et al., 2008;) - Significant improvement in patient-reported stair ambulation (Hahn et al., 2016)
Mobility	Patient has difficulty with dual tasking while walking with the prosthesis	<ul style="list-style-type: none"> - Significantly improved capacity and performance in executing a concurrent task while walking with the prosthesis (Mileusnic et al., 2017; Hahn et al., 2016; Hahn et al., 2015; Kannenberg et al., 2014; Hafner et al., 2009)

Mobility	Patient has difficulty with performing activities of daily living	<ul style="list-style-type: none"> - Significantly improved performance in the execution of various activities of daily living (Kannenberget al., 2014; Theeven et al., 2011 and 2012) - Significant improvement in PRQ Ambulation and PEQ Utility (Hahn et al., 2021) - Almost significant (p=0.056) but clinically meaningful improvement in patient-reported mobility (PLUS-M) (Davie-Smith et al., 2021)
Mobility	Patient is limited in his/her mobility Patient uses a wheelchair and a prosthesis	<ul style="list-style-type: none"> - Significant increase in over-ground walking speed of up to 25% (Hahn et al., 2021; Davie-Smith et al., 2021; Jayaraman et al., 2021; Eberly et al., 2014; Kannenberg et al., 2014; Kahle et al., 2008) - Significant improvement in distance walked in the 2-minute walk test (Davie-Smith et al., 2021) - Significant reduction in additional use of a wheelchair from 87% to 37% of subjects (Mileusnic et al., 2017) - Patients spent significantly more time active and significantly less time sitting (Kaufman et al., 2018) - About 50% of K2 patients are able to improve their overall mobility level to K3 (Hahn et al., 2021; Hahn et al., 2016; Hahn et al., 2015; Kannenberg et al., 2014; Hafner et al. 2009; Kahle et al., 2008)
Quality of life	Patient has reduced quality of life	<ul style="list-style-type: none"> - Significant improvement in health-related quality of life (Davie-Smith et al., 2021)

References

Burnfield JM, Eberly VJ, Gronely JK, Perry J, Yule WJ, Mulroy SJ. Impact of stance phase microprocessor-controlled knee prosthesis on ramp negotiation and community walking function in K2 level transfemoral amputees. *Prosthet Orthot Int* 2012;36(1):95-104.

Davie-Smith F, Carse B. Comparison of patient-reported and functional outcomes following transition from mechanical to microprocessor knee in the low-activity user with a unilateral transfemoral amputation. *Prosth Orthot Int* 2021;45(3):198-204.

Eberly VJ, Mulroy SJ, Gronley JK, Perry J, Burnfield JM. Impact of a stance phase microprocessor-controlled knee prosthesis on level walking in lower functioning individuals with transfemoral amputation. *Prosth Orthot Int* 2014;38(6):447-455.

Hahn A, Bueschges S, Prager M, Kannenberg A. The effect of microprocessor controlled exo-prosthetic knees on limited community ambulators: systematic review and meta-analysis. *Disabil Rehabil* 2021 Oct 25:1-19.

Hahn A, Lang M. Effects of mobility grade, age, and etiology on functional benefit and safety of subjects evaluated in more than 1200 C-Leg trial fittings in Germany. *J Prosthet Orthot* 2015; 27(3): 86-95.

Hafner BJ, Smith DG. Differences in function and safety between Medicare Functional Classification Level-2 and -3 transfemoral amputees and influence of prosthetic knee joint control. *J Rehabil Res Dev* 2009;46(3):417-434.

Jayaraman C, Mummidisetty CK, Albert MV, et al. Using a microprocessor knee (C-Leg) with appropriate foot transitioned individuals with dysvascular transfemoral amputations to higher performance levels: a longitudinal randomized clinical trial. *J Neuroeng Rehabil*. 2021;18(1):88.

Kahle JT, Highsmith MJ, Hubbard SL. Comparison of Non-microprocessor Knee Mechanism versus C-Leg on Prosthesis Evaluation Questionnaire, Stumbles, Falls, Walking Tests, Stair Descent, and Knee Preference; *J Rehabil Res Dev* 2008;45(1):1-14.

Kannenbergh A, Zacharias B, Pröbsting E: Benefits of microprocessor prosthetic knees to limited community ambulators: A systematic review. *J Rehabil Res Dev* 2014;51(10):1469-1495.

Kaufman KR, Bernhardt KA, Symms K. Functional assessment and satisfaction of transfemoral amputees with mobility (FASTK2): A clinical trial of microprocessor-controlled vs. non-microprocessor-controlled knees. *Clin Biomech (Bristol, Avon)* 2018 Oct;58:116-122.

Mileusnic M, Hahn A, Reiter S. Effects of a novel microprocessor-controlled knee, Kenevo, on the safety, mobility, and satisfaction of lower-activity patients with transfemoral amputation. *J Prosthet Orthot* 2017;29(4):198-205.

Lansade C, Vicaut E, Paysant J, Ménager D, Cristina MC, Braatz F, Domayer S, Pérennou D, Chiesa G. Mobility and safety with a microprocessor-controlled knee in moderately active amputees: A multi-centric randomized crossover trial. *Ann Phys Rehabil Med* 2018;61(5):278-285.

Theeven P, Hemmen B, Rings F, Meys G, Brink P, Smeets R, Seelen H. Functional added value of microprocessor-controlled knee joints in daily life performance of Medicare Functional Classification Level-2 amputees. *J Rehabil Med* 2011;43(10):906-915.

Theeven PJ, Hemmen B, Geers RP, Smeets RJ, Brink PR, Seelen HA. Influence of advanced prosthetic knee joints on perceived performance and everyday life activity of low-functional persons with a transfemoral amputation or knee disarticulation. *J Rehabil Med* 2012;44(5):454-461.

Wong CK, Rheinstein J, Stern MA. Benefits for adults with transfemoral amputation and peripheral artery disease using microprocessor compared with nonmicroprocessor prosthetic knees. *Am J Phys Med Rehabil* 2015; 94 (10): 804-810.

Ottobock North America, Reimbursement
P 800 328 4058 F 800 230 3962
US: <https://shop.ottobock.us>
CA: <https://shop.ottobock.ca>
reimbursement911@ottobock.com