

AxonSkin Silicone

As the very first prosthetic silicone glove for the Michelangelo Hand, the AxonSkin Silicone glove offers a more natural look and feel.



Description

The AxonSkin Silicone Glove is made of a multilayered silicone, therefore the features differ from the AxonSkin Natural Glove (PVC):

- Silicone is stain resistant and easy to clean
- AxonSkin Silicone does not show any change or reduced functionality in lower temperatures
- Patients will benefit from the natural look and feel. Every AxonSilicone glove is hand painted with a 2-tone coloring
- Exclusive to the Michelangelo Hand

Features/user benefits

Feature	User benefits
Hand painted details (nails, knuckles)	Very natural, aesthetically pleasing appearance, especially for user with darker skin tone
2-tone-coloring: inside of hand lighter than palm side (dark colors)	
Stain resistant, better cleaning properties	Maintains original look, easy to clean
No loss of functionality at low temperature	Reliability ensures safety in every situation

Specifications

Material	Multilayered silicone
Coating	CVD inside for better donning and doffing of glove, better interaction between system and glove
Length of sleeve	335 mm /-10mm tolerance
Circumference	197mm

Article Number

8S511=*	Men
Side	Right (R), Left (L)
Size	M
Color	2, 4, 6, 8, 10

** Coloring differs to AxonSkin Natural (PVC). Please use color swatch 646M72 (Skin Silicone Gloves) for choosing right color.

Quality for life