

1C68 Triton side flex

Your will. Your way.

The Triton family was designed for ambitious individuals who navigate varied indoor and outdoor environments and place a high value on uncompromised response and control even when performing high-impact activities.

Designed for optimal outcomes:

- Triton side flex's unique side-to-side flexibility of +/- 10° for immediate and full-ground contact while walking and standing, including on uneven surfaces
- An interconnected, linear spring design provides the level of midstance support needed for activities where responsiveness is essential
- For users with limited clearance
- Waterproof and corrosion-resistant
- Robust and maintenance-free technology
- Customizable shock absorption with two different heel wedge options
- Slim footshell option



Technical data

Suitable for	MG 3-4
Max. body weight	275 lbs / 125 kg
Sizes	22-30 cm
Footshell	Slim shape for 15 +/- 5 mm heel height (22-27 cm) Normal shape for 10 +/- 5 mm heel height (24-30 cm)
Weight without footshell*	20.6 oz / 585 g
Weight with normal footshell*	28.4 oz / 805 g
Clearance with normal footshell*	86 mm / 3 ³ / ₈ inches

* Reference size is 26 cm.



1C68 Triton side flex

Order form

Customer		Shipping address (if different from customer address)	
Customer no.	<input type="text"/>	Customer no.	<input type="text"/>
Company	<input type="text"/>	Company	<input type="text"/>
Street	<input type="text"/>	Street	<input type="text"/>
Postal code/city	<input type="text"/>	Postal code/city	<input type="text"/>
O&P professional	<input type="text"/>	Customer ref.	<input type="text"/>

Quantity	Article no.	Side	Size	- Stiffness	- P /	Color	Shape
	1C68=				- P /		
	1C68=				- P /		

Side	Size [cm]	Color	Shape
Right R	22, 23, ..., 30	Beige 4	Normal N
Left L		Light brown 15	Slim S

Selection chart *

Body weight	Sizes									
	22 cm	23 cm	24 cm	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm	
up to 121 lb (55 kg)	1	1	1	1	1	-	-	-	-	
122-165 lb (56-75 kg)	2	2	2	2	2	2	2	2	2	
166-220 lb (76-100 kg)	3	3	3	3	3	3	3	3	3	
221-275 lb (101-125 kg)	-	-	-	4	4	4	4	4	4	

Slim footshell available (15 +/- 5 mm heel height)
 Both footshell types available
 Normal footshell available (10 +/- 5 mm heel height)

* Please read the 1C68 instructions for use regarding potentially excluded combinations of configurations with Ottobock structural components.

Order example
1C68=R26-3-P/4N



The scope of delivery includes the 1C68 foot module, the 2C6 footshell including connection cover, the 2F60 heel wedge set (comprising two different degrees of hardness) and a black Spectra sock.

Recommended combinations**

TT	KD	TF
4R220* DVS & 6Y94=*-F DVS liner & 453A30/40 ProFlex Plus knee sleeve	Laminated socket & 6Y110 Skeo Sealing & 21Y21 ClickValve	7T451=1 SiOCX TF Pro; Laminated socket & 6Y110 Skeo Sealing & 21Y21 ClickValve
4R10=111 Quickchange	3C88-3* C-Leg 4 3B1*=ST Genium 3B5*=ST Genium X3	3C98-3* C-Leg 4 3B1* Genium 3B5* Genium X3
		4R57=WR/4R57=WR-ST Waterproof rotation adapter
		3F1=2 + 99B120=* Functional Cosmesis for Genium

** All components are sold separately and are available Ottobock products that are compatible with the 1C68 Triton side flex foot, which helps ensure optimal performance. Practitioners need to select components based upon individual patient criteria.

Date Place Signature