

1C10 Terion

Order Form

Company

Practitioner

Phone

Fax

Address

Email Address

Requested Due Date

PO Number

Shipping Options:

Next Day

Ground

2-Day

Other _____

Patient Information

Name _____

Age _____

Gender male female

Side of amputation left right

Weight _____ lbs

Foot size _____ cm

Stiffness selection chart

		Sizes		
		small module	medium module	large module
Body weight		22-23 cm	24-25 cm	26-28 cm
	Up to 165 lbs	Stiffness 2	Stiffness 3	Stiffness 4
	Up to 220 lbs			
Up to 275 lbs				

Slim footshell available
 Normal footshell available

Side	Size (cm)	Color	Shape
Right R	22, 23, ..., 28	Beige 4	Sizes 22 - 23cm S
Left L		Light Brown 15	Sizes 24 - 28cm N

1 Order chart for the Terion foot module including footshell

Quantity	Article no.	=	Side	Size	-	Stiffness	-	P	/	Color	Shape
	1C10	=			-		-	P	/		
	1C10	=			-		-	P	/		
	1C10	=			-		-	P	/		

Order example 1C10 = R 26 - 4 - P / 4N

2 Order chart for the Terion foot module without footshell

Quantity	Article no.	=	Side	Size	-	Stiffness	-	P	/	0
	1C10	=	N	2223	-	2	-	P	/	0
	1C10	=	N	2425	-	3	-	P	/	0
	1C10	=	N	2628	-	4	-	P	/	0

Order example 1C10 = N2628 - 4 - P / 0

3 Order chart for the footshell

Quantity	Article no.	=	Side	Size	/	Color	Shape
	2C12	=			/		
	2C12	=			/		
	2C12	=			/		

Order example 2C12 = R26 / 4 N

Additional components

Quantity	Article no.

Additional components

Quantity	Article no.