

C-Leg[®] and Kenevo[®] Selection Guide

MPK Functional Differences for K2 and K3 Patients



Features	Kenevo	C-Leg 4
Patient indication	K2; Low-K3 functional level	K3 - K4 functional level
Variable cadence	Yes; up to 1.9 mph	Yes; up to 4 mph
Stance control	Yes; varies according to chosen mode (A, B, C)	Yes; yielding resistance
Stumble recovery	Stumble Recovery Plus allows for micro-processor controlled stumble recovery activation in swing flexion and extension	Stumble Recovery Plus allows for microprocessor controlled stumble recovery activation in swing flexion and extension
Sitting function	Flexion resistance increases with increasing knee angle (progressive)	Flexion resistance can be customized for Sitting and Ramps/ Stairs descent to allow for either constant or increasing resistance during either activity
Standing function	Knee flexion blocking for support if sit-to-stand is interrupted. Intuitive stance can be activated in Mode C	Intuitive standing function is active whenever the user has the knee flexed and not moving
Stance release	Releases from stance to swing later in the gait cycle to improve stability during weight transfer	Releases from stance to swing earlier in gait cycle to facilitate weight transfer between steps
Donning function	Donning possible in the standing or seated position, as the knee can freely flex in the A, B, and B+ modes	Donning is possible from a standing or seated position
Advanced features	Wheelchair and stationary cycling for rehabilitation activities	Thousands of combinations for favorite activities with MyMode Plus
IP Rating	IP22: Designed specifically for lower mobility users. Not waterproof or corrosion-resistant, but protected from light water spray	IP67: Protected from dust, sand, dirt and temporary submersion in fresh water (up to 1m for 30min)