

4R11 Quickchange.





Amplify your possibilities.

With the 4R11 *Quickchange* adapter, pediatric and adult users can remove their distal prosthetic components from the socket by themselves when needed, in just one step. This makes dressing and undressing easier and sitting more comfortable. In addition, fitting various feet or knee-foot combinations, such as sports components, is possible.

Product features

- Four connection options allowing for flexibility of fittings – pyramid, pyramid receiver, thread and thread for modular system for children.
- The base element provides quick and easy handling for all user groups.
 - Simple one turn operation of locking mechanism
 - Arrow symbol showing open and close direction for the user
- Adapted for use with children's modular systems with addition of the 4X907=1 part
- Waterproof and corrosion resistant - protected against fresh, salt and chlorinated water.

Technical data

				
Set article number	4R11=P-1-3	4R11=R-1-3	4R11=ST-1-3	4R11=ST-2-3
Consists of (single parts)	4R11=P-1-1 4R11=R-3-2	4R11=R-1-1 4R11=R-3-2	4R11=ST-1-1 4R11=R-3-2	4R11=ST-2-1 4R11=R-3-2
Field of application	Adults/children: TT, KD, TF, HD (only adults)			
K-Level/ Mobility grade	K1 - K4 / MG1 - MG4			
Max. body weight	330 lbs / 150 kg			
Max. body weight for sports	275 lbs / 125 kg			
Build height	1.73" (44 mm)	2.17" (55 mm)	1.61" (41 mm)	1.61" (41 mm) / 1.69" (43 mm) with 4X907=1
Weight	13.05 oz (370 g)	13.05 oz (370 g)	11.64 oz (330 g)	11.46 oz (325 g)
Connection	Pyramid/ pyramid receiver	Pyramid receiver/ pyramid receiver	Thread/ pyramid receiver	Thread for modular system for children/ pyramid receiver



4R11 Quickchange. Order form.

Customer

Customer no.

Company

Street

Zip/Postal code City

O&P professional

Shipping address (if different from customer address)

Customer no.





Company

Street







Zip/Postal code City

Customer

Set (consisting of 1 connecting element and 1 base element)

Quantity	Article no.	Description	
	4R11=P-1-3	Quickchange set with pyramid	
	4R11=R-1-3	Quickchange set with pyramid receiver	
	4R11=ST-1-3	Quickchange set with thread	
	4R11=ST-2-3	Quickchange set with thread for modular system for children	

Single parts and accessories

Quantity	Article no.	Description	
	4R11=R-3-2	Base element with pyramid receiver	
	4R11=P-1-1	Connecting element with pyramid	
	4R11=R-1-1	Connecting element with pyramid receiver	
	4R11=ST-1-1	Connecting element with thread	
	4R11=ST-2-1	Connecting element with thread for modular system for children	
	4X907=1	Reduction pyramid receiver for modular system for children	

Application examples

Easy removal of the prosthetic components below the socket for comfortable sitting/
to change clothes more easily (1 set)



Change from everyday to secondary prosthesis (1 set + 1 additional base element)



Change between two sockets (1 set + 1 additional connecting element)

