

6Y112 Skeo Sealing TT.

The 6Y112 *Skeo Sealing TT* liner is the latest addition to our Skeo Sealing Liner family, offering transtibial amputees a high level of comfort and reliability. The sealing ring provides a consistent seal against the inner socket wall for secure suspension. Its silky-smooth surface with no textile cover simplifies cleaning as well as the donning and doffing process.

Designed for optimal outcomes

- One sealing height provides a higher vacuum area for higher adhesion as well as security and peace-of-mind for the patient
- Integrated sealing ring promotes durability throughout the lifetime of the liner
- The inside of the liner is textured for a pleasant feeling on the skin
- Silky-smooth outer coating for easier donning and doffing
- Orientation aid on the inside of the liner for an optimal liner position on the residual limb
- An integrated matrix for reduced elongation and good tissue management
- The *Skeo Sealing TT* can be combined with the NEW 21Y2 *QuickValve* or the *Harmony* system plus 21Y2 *QuickValve*.
- Waterproof: Protected against fresh, salt and chlorinated water



6Y112

Skeo Sealing TT.


| Article number | Circumference | Ring height |
|----------------|---------------|-------------|
| 6Y112=160X80 | 160 mm | 80 mm |
| 6Y112=180X80 | 180 mm | 80 mm |
| 6Y112=200X80 | 200 mm | 80 mm |
| 6Y112=210X80 | 210 mm | 80 mm |
| 6Y112=220X80 | 220 mm | 80 mm |
| 6Y112=235X80 | 235 mm | 80 mm |
| 6Y112=250X80 | 250 mm | 80 mm |
| 6Y112=265X80 | 265 mm | 80 mm |
| 6Y112=280X80 | 280 mm | 80 mm |
| 6Y112=300X80 | 300 mm | 80 mm |
| 6Y112=320X80 | 320 mm | 80 mm |
| 6Y112=340X90 | 340 mm | 90 mm |
| 6Y112=360X90 | 360 mm | 90 mm |
| 6Y112=380X90 | 380 mm | 90 mm |
| 6Y112=400X90 | 400 mm | 90 mm |



The ring



The material



The orientation aid