

21Y2 QuickValve.

The 21Y2 *QuickValve* is an expulsion valve and can be opened manually. It is intended for fabricating vacuum socket systems, for example in combination with a Skeo Sealing liner or all common TT liners without a distal connection in combination with a sealing sleeve.

Designed for optimal outcomes

- Wide range of applications in combination with a liner (TT, KD, TF)
- The valve can be easily opened by turning the valve cap with recessed grips
- When closing, a noticeable signal indicates that the valve is closed
- The orientation rectangles indicate where the valve cover is placed on the valve base
- The M10 thread, allows for easy installation into the socket
- Easy to clean
- Waterproof protected against fresh, salt and chlorinated water



Article number	21Y2=M10
Mobility grade	1 – 4
Amputation level	TT, TF, KD
Max. body weight	No limit
Total height of the valve body	0.66" (17 mm)
External height (overall height outside the socket) of the valve body	0.43" (11 mm)
Outer diameter of the valve cap	1.10" (28 mm)

