

ottobock.

MagnoFlex Lock

The shuttle lock with
magnetic pin guide



Quality for life

Information for users

Easy to put on – the MagnoFlex Lock with flexible pin

A shuttle lock allows a firm connection between the residual limb and the prosthesis. A connection you can establish much easier in the future with the MagnoFlex Lock, the first to combine a flexible pin and a shuttle lock with integrated magnet. The result: a straightforward pin guide, making it easier and simpler than ever before to put on your prosthesis.





• **MagnoFlex Lock**
with flexible pin



A strong connection – with many advantages

An innovation you will appreciate: the straightforward pin guide

The MagnoFlex Lock consists of a flexible pin and a housing made of high-performance synthetic material, as well as ratchet unit, which set themselves apart with a very low weight and high durability.

A magnet which entirely surrounds the guide opening for the socket attachment block is integrated into the housing. Thanks to the magnet, the pin is automatically attracted and the connection between the residual limb and the prosthesis can be made quickly and easily. The pin guide is supported by the innovative pin: because it is flexible, the pin, attracted by the magnet, simply slides along the conical sidewalls of the socket attachment to the correct position and locks in place.

That something extra – included

With the MagnoFlex Lock, you can have your fitting converted quickly and easily – on a trial or permanent basis – from a mechanical to a suction fitting in order to also benefit from the advantages of suction technology. For test purposes, the flexible pin is unscrewed from the liner and the opening in the socket attachment is sealed with the valve included in the scope of delivery – finished.

Your benefits at a glance

- Thanks to the straightforward pin guide, the liner engages in the prosthesis quickly and easily – without arduous, time-consuming searching.
- The use of high-performance synthetic material reduces the weight of the prosthesis while offering great durability.
- With the enclosed valve, the MagnoFlex Lock can be converted to a suction fitting quickly and easily by your practitioner – temporarily or permanently.

With compliments from



Otto Bock HealthCare GmbH
Max-Näder-Straße 15 · 37115 Duderstadt/Germany
T +49 5527 848-1706 · F +49 5527 72330
export@ottobock.de · www.ottobock.com