

ottobock.



4X880*
Genium Protector

Quality for life

4X880* Genium Protector

Protection – Volume – Attractive Design



Visualised knee rotation point
Orientation aid for prosthetic alignment

Easy to use closure technology
Allows user to put on and take off the main part of the Protector

Natural calf shape
For a less conspicuous appearance underneath clothing

Recess for inductive charger
Allows for inductive charging even through clothing

Thick wall
Specialized protection for the Genium charging unit

Anatomically-shaped knee section
Facilitates kneeling

Smooth surface
Makes it easier to dress and undress

Tough main component
Protects the Genium prosthesis system from jolts, the environment and normal wear and tear

Area for length adjustment
The main Protector component can be shortened by qualified personnel

Anti-slip fabric
Helps keep socks from slipping off

Innovative foot cuff
Hybrid construction of plastic and fabric offers flexibility and maintains shape

The Genium Protector basically consists of the main protector component with an innovative footshell and protects the Genium prosthesis system against impacts, environmental influences and wear. In addition to protection, the Genium Protector also provides good volume equalisation under clothing. Furthermore it offers a variety of daily functions like inductive charging option of the Genium and a specially designed knee area to facilitate kneeling.

Technical data

Reference number	4X880*
Material	Plastic/Textile
Weight (Protector)	450 g
Weight (foot cuff)	~ 60 g
Sizes (incl. foot cuff)	S, M, L
Colour	Champaign/black

Compatible prosthetic feet

Prosthetic feet*		Foot cuff		
Product Name	Reference number	Size S	Size M	Size L
Triton	1C60	24–25	26–28	29–30
Triton Vertical Shock	1C61	24–25	26–28	29–30
Triton Low Profile	1C63	24–25	26–28	29–30
Triton Heavy Duty	1C64	24–25	26–28	29–30
Trias	1C30	24–25	26–28	29–30
C-Walk	1C40	24–25	26–28	29–30
Dynamic Motion	1D35	23–25	26–28	29–30
Adjust	1M10	23–25	26–28	29–30
Axition	1E65	23–25	26–28	29–31
Lo Rider	1E57	24–25	26–28	29–31

Order example

Reference number	=	Consisting foot cuff size
4X880	=	M

*The foot cuff is not compatible with a slim footshell.