








































# Foot Family Portfolio.

ottobock.

Product Image	K-Level	Material Number	Foot Name	Foot Description	Max Pt Weight	Clearance - Inches (size 26)	L-Codes
	K3	1A1-2	Empower®	Microprocessor Feet	287 lbs	8 3/4" (222 mm)	L5973 + L5969
	K3	1B1-2	Meridium®	Microprocessor Feet	275 lbs	6 1/4" (160 mm)	L5973 
	K3	1C50	Taleo	Moderate-High Activity	330 lbs	5 7/8" (150 mm)	L5981 
	K3	1C51	Taleo Vertical Shock	Moderate-High Activity	330 lbs	7 1/4" (185 mm)	L5987 + L5986 
	K3	1C52	Taleo Harmony	Moderate-High Activity	330 lbs	7 1/4" (185 mm)	L5987 + L5986 + L5781 
	K3	1C53	Taleo Low Profile	Moderate-High Activity	330 lbs	2 1/2" (65 mm)	L5981
	K3	1C56	Taleo Adjust	Moderate-High Activity	253 lbs	4 5/8" (118 mm)	L5981 + L5990 
	K3	1C58	Taleo Side Flex	Moderate-High Activity	287 lbs	3 13/16" (97 mm)	L5981 + L5986 
	K3-K4	1C60	Triton®	Moderate-High Activity	330 lbs	6 1/4" (158 mm)	L5981 
	K3-K4	1C61	Triton® Vertical Shock	Moderate-High Activity	330 lbs	8 1/16" (205 mm)	L5986 + L5987 
	K3-K4	1C62	Triton® Harmony	Moderate-High Activity	330 lbs	8 1/16" (205 mm)	L5986 + L5987 + L5781 
	K3-K4	1C63	Triton® Low Profile	Moderate-High Activity	330 lbs	2 7/8" (77 mm)	L5981
	K3-K4	1C64	Triton® HD	Moderate-High Activity	330 lbs	5 7/8" (149 mm)	L5981 
	K3-K4	1C68	Triton® Side Flex	Moderate-High Activity	275 lbs	3 3/8" (86 mm)	L5981 + L5986 
	K3-K4	1E57	Lo Rider™	Moderate-High Activity	300 lbs	1 3/8" (36 mm) (with pyramid)	L5981
	K3-K4	1C20	ProSymes™	Specialty	275 lbs	1 3/4" (43 mm)	L5981
	K3	1C30	Trias®	Low-Moderate Activity	275 lbs	4 3/8" (111 mm)	L5981 + L5986
	K3-K4	F23	Maverick Vertical Shock	Moderate-High Activity	324 lbs	8 1/2" (216 mm)	L5987 + L5986 
	K3-K4	F21	Maverick Xtreme AT	Moderate-High Activity	365 lbs	6 11/16" (170 mm)	L5981
	K3	F22	Maverick Comfort AT	Low-Moderate Activity	365 lbs	4 5/8" (117 mm)	L5981
	K3	1D35	Dynamic Motion™	Moderate-High Activity	220 lbs	3 3/8" (87 mm)	L5979
	K2	VS4	Kintrol	Lower Activity	300 lbs	4 7/16" (112 mm)	L5972 + L5968
	K1-K2	1C11	Terion™ K2	Lower Activity	385 lbs	3" (77 mm)	L5972 + L5986
	K2	VS5	Restore	Lower Activity	300 lbs	3 1/2" - 3 7/8" (89 - 99 mm)	L5972 + L5986
	K1-K2	1D10/ 1D11	Dynamic Foot™ with Toes	Lower Activity	330 lbs	3 1/2" (90 mm) (without adapter)	L5972
	K1-K2	1S101/102/103	SACH+™ Feet	Lower Activity	275 lbs	3 3/8" (85 mm)	5971 (Replacement) Included in base code
	K1-K3	1WR95	Foot for water resistant prosthesis	Specialty	330lbs	3 3/8" (85 mm)	L5972